

Kinematics Lab – Day 1: Measurement

Utley, A. (2019), Chapter 3: *Measurement and assessment in motor control*: In: *Motor control, learning and development: Instant notes* (Second edition). Routledge.

Ansuini, C., Giosa, L., Turella, L., Altoè, G., & Castiello, U. (2008). An object for an action, the same object for other actions: Effects on hand shaping. *Experimental Brain Research*, 185(1), 111–119.
<https://doi.org/10.1007/s00221-007-1136-4>

Today's topics

- General on motion tracking
- Transforming Ansuini et al. into our demo study
- Data collection of one very small and one proper demo

Today's schedule

8:15 – 8:45	Lecture: Introduction to the lab, equipment, and motion capture	
9:00 – 13.00	Motion lab with Cordula & Silke; Group-work with Laura	
	9:00	Group 1
	9:30	Group 2
	10:00	Group 3
	10:30	Group 4
	11:00	Group 5
	11:30	Group 6
	12:00	Group 7
	12:30	Group 8
	13:00	Group 9 (?)
13:30 – 14:00	Classroom: Wrap-up	

Best not to even bring any rings, watches, wristbands etc. –
Nothing shiny!

Also: No food or drinks

Behave carefully and considerate, esp. given the small space

Read the exercise instructions before coming to the lab

The value of general (marker-based or marker-less) motion tracking

Motion tracking refers to all kinds of techniques to record dynamic behavior in space.

Why use motion tracking?

- Measure very small cues hidden in actions, e.g. about the underlying intention
- Quantify various performance parameters that are hidden in e.g. RT measures
- Precise analyses of movements, e.g. in sports sciences, training, and rehabilitation research
- Measure task-specific deviations from 'typical' performance
- Allow online control of experimental procedure based on action performance
- Use more realistic action tasks (in contrast to virtual movements etc.)
- Control that no movements happened, e.g. in infant research
- Detect early neuropsychological problems, e.g. in Autism

Deciding on measurement technique(s)

Decisions to make:

- Which parts of the body (single limbs, a finger, an arm, the whole body)
- What spatial coverage (one or multiple sensors, placement of sensors)
- What technology (optical, magnetic, gravity-based), what gets measured (position, orientation, acceleration ...)
- What temporal and spatial resolution (high, low)
- How prone to occlusions and other noise
- How to ensure participants' comfort and safety (wearable, not disturb natural body movements and task performance)
- How portable (stationary systems, more flexible systems, or simple devices like mobile phones)
- Many more practical aspects: price range and availability, software availability, ease of use, room requirements etc.

Tracking systems

Specific movement tracking devices

OptoTrac, Vicon, Qualisys

Polhemus

Video camera based

Sensor types

Sensors (w/o cables)

Body suits

Sensor gloves

Commercial performance trackers

Wearable sensors (google watch etc.)

Objects with built-in sensors (mobile phones)

Gaming technology (often sensor-less)

Kinect, Wii

Virtual Reality

Motion data usually require preprocessing; afterwards they can be treated like any other data.

Tracking systems: Magnetic sensors

How does it work?

Source creates light magnetic environment

Sensors detect distortions due to movement

Advantages

Not so sensitive to occlusions

Flexible setup, portable


Disadvantages

Requires shielded environment

Difficult to combine with other methods, e.g. EEG



Example: Polhemus G4 system



Tracking systems: Optical sensors

How does it work?

Active or passive sensors are picked up with a system of cameras

Common advantages

Stable signal

Large coverage

Common disadvantages

Prone to occlusions, requires many cameras

Normally requires fixed setup

Other systems, e.g. Microsoft Kinect

Calibration procedure to detect joints → extracts skeleton based on movement contingencies: no sensors required



Example: Vicon system

Tracking systems: Accelerometers

How does it work?

Inertia-based, measures force against gravity

Advantages

Position independent

Often easy to set up

Disadvantages

Does not inform about absolute spatial position

Examples

Smart watches

Mobile phones

Specialized equipment, e.g. Swimtag



Example: Swimtag

The IMC Motion Lab



MoTEC lab

(Motion Tracking for Embodied Cognition)

Room 1483-433

Access with keycard

Funding by:

It-vest
samarbejdende universiteter

Optical motion tracking: Qualisys Miquis

www.qualisys.com

Optical marker-based tracking

8 Miquis M3 infrared cameras

1 synchronized Miquis video camera

Passive, reflective markers


Constant sample rate of 100 Hz

Tracking of many markers possible

QUALISYS



6 degrees of freedom



19 mm



16 mm



14 mm



12.5 mm




9.5 mm



8 mm



6.5 mm




How difficult is it to measure movement at IMC Motion Lab?

Easy 😊

- Full setup with computer, software, cameras, markers etc.
- Quite intuitive software package: Qualysis Track Manager (QTM)
- Ready-made arrangement of cameras in the IMC Motion Lab
- Intuitive measures

Challenging 😞

- Requires thorough planning, testing and learning new techniques
- Risk of data loss or marker switching
- Analysis and interpretation of movement not always trivial
- Requires really good experiment designs



Navigating today's lab visit

Qualisys Track Manager (QTM)

Matlab for streaming and experiment control

Automatic Identification of Markers (AIM) models

Many preparations done already

Transforming Ansuini et al. into our demo study

“The aim of the present study was to investigate whether the prior-to-contact grasping kinematics of the same object vary as a result of different goals of the person grasping it.” (Ansuini et al., 2008, p. 111)

Conditions

- (1) **Grasp** it without performing any subsequent action
- (2) Lift and **throw** it
- (3) **Pour** the water into a container
- (4) **Place** it accurately on a target area
- (5) **Pass** it to another person

Our variables:

MovementOnset, reachDuration, peakReachVelocity, maxReachHeight, maxReachAperture, graspHeight

Today's schedule

8:15 – 8:45	Lecture: Introduction to the lab, equipment, and motion capture		
9:00 – 13.00	Motion lab with Cordula & Sille; Group-work with Laura		
	9:00	Group 1	
	9:30	Group 2	
	10:00	Group 3	Best not to even bring any rings, watches, wristbands etc. – Nothing shiny!
	10:30	Group 4	
	11:00	Group 5	
	11:30	Group 6	Also: No food or drinks
	12:00	Group 7	
	12:30	Group 8	Behave carefully and considerate, esp. given the small space
	13:00	Group 9 (?)	
13:30 – 14:00	Classroom: Wrap-up		